

SeraCare Life Sciences
Human Papillomavirus Genotype Performance Panel
PSH301-Data Sheet

Table 1. Typical Data for Genotype Performance Panel, PSH301, on Various Test Platforms

Test Methods	Vendor	HPV Subtype	Result
hc2 Hybrid Capture® 2 High-Risk HPV DNA Test®	Qiagen	HPV16	High Risk Positive
		HPV31	High Risk Positive
		HPV18	High Risk Positive
		Negative Control	--
hc2 Hybrid Capture® 2 HPV DNA Test® (Low-Risk HPV Probe)	Qiagen	HPV11	Low Risk Positive
		Negative Control	--
Linear Array® HPV Genotyping Test	Roche Diagnostics	HPV16	16
		HPV11	11
		HPV31	31
		HPV18	18
		Negative Control	HPV DNA Not Detected
Cervista™ HPV HR Test	Hologic, Inc.	HPV16	Positive
		HPV31	Positive
		HPV18	Positive
		Negative Control	Negative
Cervista™ HPV 16/18 Test	Hologic, Inc.	HPV16	HPV16
		HPV18	HPV18
HPV HR ASR	Hologic, Inc.	HPV16	A9 Oligo Positive
		HPV31	A9 Oligo Positive
		HPV18	A7 Oligo Positive
		Negative Control	--

For additional information, please contact SeraCare Life Sciences Technical Support

Phone (508)244-6400 or US Toll Free (800) 676-1881

Visit our website at www.seracare.com or email us at info@seracare.com



25 Birch Street
Milford, MA 01757 USA

April 2009
12090-02